Assigned Values	REF W.MALB.CTL	
	LOT D121	
ANALYTE		mg / l
Micro Albumin (urine)	mALB	31 - 46

BIOMEDICINA - PATOLOGIA CLINICA - INFORMATICA

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## INTENDED USE

This product line (Controls) is for use in Quality Control procedures in the determination of:

- **Micro Albúmin (mALB)**, in human urine carried out with NSC's kits.

## METHOD

The method is based on the Immunoprecipitation reaction in liquid phase (turbidimetry) with specific antibodies. For further information, see the instructions leaflet included in the kits.

## CONTENTS - COMPONENTS

- *Micro Albumin (urine) - Low Control (W.MALB.CTL)* 6 vials of 1 ml

**Form:** liquid, ready for use.

**Stability:** 12 months at +2°C ... +8°C.

## INFORMATION - PREPARATION

**Origin:** Human serum.

Due to lack of provisions from the Scientific Committees, our "Internal Standardisation" has been applied.

### Concentration

See table on the back.

## PRECAUTIONS AND HAZARDS

- **Preservative: Sodium Azide.**

- **Handling:** Sodium Azide is toxic.

All the reagents contain Sodium Azide <0,1% (w/v). Even if Sodium Azide is not dangerous at this concentration, make sure to take the necessary precautions to prevent accidental swallowing.

- **Waste:** Sodium Azide can react with lead and copper plumbing to form highly explosive metal azides.

- **Calibrators and Controls – Serum and Plasma of human origin.**

- **Calibrators and Controls Supplier Statement:**

Serum or plasma of human origin has been tested by an "FDA – USA" approved method and found negative for the presence of Hepatitis B surface antigen (HBsAg), HCV antibodies, and HIV-1, 2 antibodies.

- **Handling and waste:**

Because no test method can completely exclude the presence of infectious agents, all materials of human origin should be handled as recommended for any potentially infectious human serum or blood specimen following the precaution recommended for biohazard materials.

## WARNINGS

- All reagents are intended only for "*In Vitro* research".
- Reagents have passed quality control checks and are expected to react as described in this Instructions. If controls do not yield the preveded reaction, all the reagents would not be considered reliable.
- Do not mix reagents from different lot kits.
- Calibrator dilutions are stable for a few hours.

## TECHNICAL INFORMATION

Contact Customer Support Service or visit [www.newscientific.com](http://www.newscientific.com).

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