NEW SCIENTIFIC COMPANY España, S.r.l.

Valencia, 558 - 08026 Barcelona

Teléfono: (+34) 93 244 82 94

Facsimil: (+34) 93 244 82 95

e-mail: nsce@newscientific.com

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1. **Solomon A**. Bence-Jones Proteins and Light Chains of Immunoglobulins (First of Two Parts). *N Eng J Med*, Vol. 294 No. 1, Jan. 1976
2. **Solomon A**. Bence-Jones Proteins and Light Chains of Immunoglobulins (Second of Two Parts). *N Eng J Med*, Vol. 294 No. 2, Jan. 1976
3. **Pezzoli A, Pascali E**. The Clinical Significance of Pure Bence Jones Proteinuria at Low Concentration. *A.J.C.P.*, Apr. 1989
4. **Solomon A**. Bence-Jones Proteins: Malignant or Benign?. *N Eng J Med*, Vol. 306 No. 10, Mar. 1982
5. **Whicher JT, Calvin J, Riches P, Warren C**. The Laboratory Investigation of Paraproteinemia. *Ann. Clin. Biochem*, 24: 119-132, 1987
6. **Stevens FJ, Solomon A, Schiffer M**. Bence-Jones Proteins: A powerful Tool for the Fundamental Study of Protein Chemistry and Pathophysiology. *Biochemistry*, Vol. 30 No. 28, 1991
7. **Matsuda K, Hiratsuka N, Koyama T, ...**. Sensitive Method for Detection and Semiquantitation of Bence Jones Protein by Cellulose Acetate Membrane Electrophoresis Using Colloidal Silver Staining. *Clinical Chemistry*, 47, no. 4, Apr. 2001
8. **Kyle R, Greipp P**. "Idiopathic" Bence Jones Proteinuria: Long-Term Follow-up in Seven Patients. *N Eng J Med*, Vol. 306 No. 10, Mar. 1982
9. **Tillyer CR**. Clinical Applications of Immunoglobulin Free Light Chain Estimations. *Int J Clin Lab Res*, 23(1), 1993
10. **Okada K, Oguchi N, Shinohara K, ...**. BUN, Bence Jones protein, and chromosomal aberrations predict survival in multiple myeloma. *Rinsho Ketsueki*, 38(12), Dec. 1997
11. **Vegh Z, Otto S, Eckhardt S**. Monoclonal Free Light Chains in Urine and their significance in clinical diagnostics: are they really tumor markers?. *J Clin Lab Anal*, 4(6), 1990
12. **Singhal S, Mehta J, Desikan R, et al**. Antitumor Activity of Thalidomide in Refractory Multiple Myeloma. *N Engl J Med*, 341(21), 1999
13. **Shinohara H, Matsuura K**. Does catalytic activity of Bence-Jones proteins contribute to the pathogenesis of Multiple Myeloma?. *Appl Biochem Biotechnol*, 83(1-3), Jan.-Mar. 2000
14. **Spasov E, Goranova V**. Prognostic assessment of the Durie and Salmon staging system in patients with multiple myeloma. *Folia Med (Plovdiv)*, 40(3B suppl 3), 1998
15. **Abe M, Goto T, Kosaka M, ...**. Differences in kappa to lambda (kappa:lambda) ratios of serum and urinary free light chains. *Clin Exp Immunol*, 111(2), Feb. 1998
16. **Uchida M, Kamata K, Okubo M**. Renal dysfunction in multiple myeloma. *Intern Med*, 34(5), May 1995
17. **Goranova V, Spasov E**. Prognostic significance of the immunological variant in patients with multiple myeloma. *Folia Med (Plovdiv)*, 41(1), 1999
18. **Weber DM, Dimopoulos MA, Mouloupoulos LA, ...**. Prognostic features of asymptomatic multiple myeloma. *Br J Haematol*, 97(4), Jun. 1997
19. **Kadokura M, Tanio N, Nonaka M, ...**. A surgical case of solitary plasmacytoma of rib origin with biclonal gammopathy. *Jpn J Clin Oncol*, 30(4), Apr. 2000
20. **Knudsen LM, Hjorth M, Hippe E**. Renal failure in multiple myeloma: reversibility and impact on the prognosis. Nordic Myeloma Study Group. *Eur J Haematol*, 65(3), Sep. 2000
21. **van Zaanen HC, Diderich PP, Pegels JG, ...**. Renal insufficiency due to light chain multiple myeloma. *Ned Tijdschr Geneesk*, 144(45), Nov. 2000
22. **Vivaldi P, Comotti C, Pedrazzoli M**. The kidney in multiple myeloma. The physiopathological and clinical aspects. *Recenti Prog Med*, 85(2):123-33, Feb. 1994
23. **Tillyer CR, Iqbal J, Raymond J, ...**. Immunoturbidimetric assay for estimating free light chains of immunoglobulins in urine and serum. *J Clin Pathol*, 44(6), Jun. 1991
24. **Pascali E**. Clinical applications of immunoglobulin free light chain analysis. *Int J Clin Lab Res*, 24(2), 1994
25. **Ivanyi B**. Development of chronic renal failure in patients with multiple myeloma. *Arch Pathol Lab med*, 117(8), Aug. 1993
26. **Irish AB, Winearls CG, Littlewood T**. Presentation and survival of patients with severe renal failure and myeloma. *QJM*, 90(12), Dec. 1997
27. **Sakhujia V, Jha V, Varma S, ...**. Renal involvement in multiple myeloma: a 10-year study. *Ren Fail*, 22(4), 2000
28. **Stevens FJ, Kisilevsky R**. Immunoglobulin light chains, glycosaminoglycans, and amyloid. *Cell Mol Life Sci*, 57(3), Mar 2000
29. **Messiaen T, Deret S, Mougenot B, ...**. Adult Fanconi syndrome secondary to light chain gammopathy. Clinicopathologic heterogeneity and unusual features in 11 patients. *Medicine (Baltimore)*, 79(3), May 2000
30. **Pote A, Zwizinski C, Simon EE, ...**. Cytotoxicity of myeloma light chains in cultured human kidney proximal tubule cells. *Am J Kidney Dis*, 36(4), Oct. 2000
31. **Doi K, Teramoto S, Hosoi T, ...**. Renal tubular acidosis type II secondary to gamma-light chain excretion in an elderly patient with multiple myeloma. *Nippon Ronen Igakkai Zasshi*, 35(6), Jun. 1998
32. **Cohen G, Haag-Weber M, Horl WH**. Immune dysfunction in uremia. *Kidney Int Suppl*, 62:S79-82, Nov. 1997
33. **Deret S, Denoroy L, Lamarine M, ...**. Kappa light chain-associated Fanconi's syndrome: molecular analysis of monoclonal immunoglobulin light chains from patients with and without intracellular crystals. *Protein Eng*, 12(4), Apr. 1999
34. **Picken MM, Shen S**. Immunoglobulin light chains and the kidney: an overview. *Ultrastruct Pathol*, 18(1-2), Apr 1994
35. **Keren DF**. Procedures for the evaluation of monoclonal immunoglobulins. *Arch Pathol Lab Med*, 123(2), Feb. 1999
36. **Qiu W, Wu G**. Immunological diagnosis of light chain disease: analysis of 11 cases found in Fujian Province. *Chung Kuo I Hsueh Yuan Hsueh Pao*, 15(4), Aug. 1993
37. **Bradwell A, Carr-Smith H, Mead G, ...**. Highly Sensitive, Automated Immunoassay for Immunoglobulin Free Light Chains in Serum and Urine. *Clinical Chemistry*, 47:4, Apr. 2001
38. **Fang L**. Light Chain Nephropathy. *Kidney International*, vol. 27, 1985
39. **Tillyer CR**. The estimation of free light chains of immunoglobulins in biological fluids. *Int J Clin Lab Res*, 22:152-158, 1992
40. **Molby L, Hansen HH, Jensen EL**. Development and treatment of renal insufficiency in multiple myeloma. *Ugeskr Laeger*, 156(30):4343-7, Jul 1994

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9. Clinical applications of immunoglobulin free light chain estimations

Int J Clin Lab Res 1993;23(1):25-9 (ISSN: 0940-5437)

Tillyer CR

Department of Chemical Pathology, Royal Marsden Hospital, London, UK.

The relevance of free light chain assays to diagnosis, staging, treatment and prognosis assessment in B-cell disorders (including myeloma, chronic lymphocytic leukaemia, and lymphoma), multiple sclerosis and diabetes is discussed and their actual and potential use is examined.

10. BUN, Bence Jones protein, and chromosomal aberrations predict survival in multiple myeloma

Rinsho Ketsueki 1997 Dec;38(12):1254-62 (ISSN: 0485-1439)

Okada K; Oguchi N; Shinohara K; Tamura N; Ishii K; Noguchi Y; Hayashi S; Yamamoto H; Takeichi M; Fujimoto H; Shirota T; Hayashi T

Third Department of Internal Medicine, Tokyo Medical College.

The prognostic significance of some clinical features in 41 patients with multiple myeloma (including 2 patients with plasma cell leukemia) from our institution was analyzed. Out of 14 variables isolated from the univariate analysis ($P < 0.05$), only three (BUN, Bence Jones protein, and chromosomal aberrations) were significant in the multivariate model ($P < 0.05$). Derived from these three variables, three subpopulations of patients were identified. The first group included 20 patients with a low risk of death and their median survival has not been reached. In particular, no one died during the first 60 months in this group. The second group also included 14 patients with an intermediate risk of death and a median survival of 49.2 months. The third group comprised seven patients with a high risk of death during 24 months after diagnosis and a median survival of 31.6 months ($P < 0.0001$). Finally, Durie & Salmon's myeloma staging system was demonstrated in the present series, and it showed prognostic validity for each stage ($P < 0.0222$). Compared with Durie & Salmon's staging system, our prognostic model for multiple myeloma was more useful to predict prognosis when applied to the present series.

11. Monoclonal free light chains in urine and their significance in clinical diagnostics: are they really tumor markers?

J Clin Lab Anal 1990;4(6):443-8 (ISSN: 0887-8013)

Vegh Z; Otto S; Eckhardt S

National Institute of Oncology, Budapest, Hungary.

Bence Jones proteins (monoclonal free light chains of immunoglobulins) are the earliest known biological markers of malignant cell dyscrasia; Bence Jones proteinuria is also present in many types of B cell-related neoplasms. Sometimes, it may also occur in Hodgkin's disease. In some cases, benign monoclonal gammopathy was found to be associated nontumorous diseases as well. The type of monoclonal light chain, the degree of polymerization, and the isoelectric point of the molecule may affect the course of the disease. Urine samples from 637 patients with true or suspected lymphoproliferative diseases were investigated over a 2-yr period by different immunochemical methods. Bence Jones proteinuria was identified in 71 cases by isoelectric focusing combined with immunofixation, while the pathological protein was detected only in 63 cases by conventional methods. Bence Jones proteins can be detected by this new method at a level below the sensitivity of conventional procedures. Bence Jones proteins in the urine may signal a malignant tumor or malignant transformation of an earlier disease. The early detection of monoclonal immunoglobulin light chains in the urine may be important in clinical diagnosis, therapy, and follow-up.

12. Antitumor activity of thalidomide in refractory multiple myeloma

N Engl J Med 1999 Nov 18;341(21):1565-71 (ISSN: 0028-4793)

Singhal S; Mehta J; Desikan R; Ayers D; Roberson P; Eddlemon P; Munshi N; Anaissie E; Wilson C; Dhodapkar M; Zeddis J; Barlogie B

Myeloma and Lymphoma Program, South Carolina Cancer Center, University of South Carolina, Columbia, USA.

BACKGROUND: Patients with myeloma who relapse after high-dose chemotherapy have few therapeutic options. Since increased bone marrow vascularity imparts a poor prognosis in myeloma, we evaluated the efficacy of thalidomide, which has antiangiogenic properties, in patients with refractory disease.

METHODS: Eighty-four previously treated patients with refractory myeloma (76 with a relapse after high-dose chemotherapy) received oral thalidomide as a single agent for a median of 80 days (range, 2 to 465). The starting dose was 200 mg daily, and the dose was increased by 200 mg every two weeks until it reached 800 mg per day. Response was assessed on the basis of a reduction of the myeloma protein in serum or Bence Jones protein in urine that lasted for at least six weeks.

RESULTS: The serum or urine levels of paraprotein were reduced by at least 90 percent in eight patients (two had a complete remission), at least 75 percent in six patients, at least 50 percent in seven patients, and at least 25 percent in six patients, for a total rate of response of 32 percent. Reductions in the paraprotein levels were apparent within two months in 78 percent of the patients with a response and were associated with decreased numbers of plasma cells in bone marrow and increased hemoglobin levels. The microvascular density of bone marrow did not change significantly in patients with a response. At least one third of the patients had mild or moderate constipation, weakness or fatigue, or somnolence. More severe adverse effects were infrequent (occurring in less than 10 percent of patients), and hematologic effects were rare. As of the most recent follow-up, 36 patients had died (30 with no response and 6 with a response). After 12 months of follow-up, Kaplan-Meier estimates of the mean (+/-SE) rates of event-free survival and overall survival for all patients were 22+/-5 percent and 58+/-5 percent, respectively.

CONCLUSIONS: Thalidomide is active against advanced myeloma. It can induce marked and durable responses in some patients with multiple myeloma, including those who relapse after high-dose chemotherapy.

21. Renal insufficiency due to light chain multiple myeloma

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In 3 female patients, aged 65, 83 and 76 years, with severe renal failure, light chain multiple myeloma was diagnosed, following a substantial delay on the part of the doctors concerned. Either the diagnosis had not suspected or the serum proteins had been misinterpreted. After a while, the first two patients declined further treatment with chemotherapy and haemodialysis, and subsequently died. The third patient attained a creatinine clearance of 20 ml/min and was subsequently treated for the multiple myeloma in the outpatients department. The absence of a paraprotein peak in the serum does not exclude the possibility of a multiple myeloma. In the case of light chain disease, the gammaglobulin region is, in fact, often empty. Treatment of multiple myeloma consists of a rapid rehydration and forced diuresis; the usefulness of plasmapheresis has not been demonstrated.

22. The kidney in multiple myeloma. The physiopathological and clinical aspects

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The renal concern in a multiple myeloma (MM) case has a frequency of 50% and causes a worsening of the disease with a survival average of about 12 months. The monoclonal light free chains (CLL) produced in excess by the plasmacytes are present in the urine as proteinuria of Bence Jones (PBJ) in 60-70% of patients affected by MM. They represent the major pathogenetic factor of the nephropathy in course of MM as they can deposit in shape of intratubular "casts" in the myeloma casts nephropathy (MCN). In some worse cases, dehydration or hypercalcaemia can cause an irreversible acute renal insufficiency (RI). It is therefore important in a patient affected with MM with PBJ to prevent, locate and opportunely treat these situations which worsen the nephropathy. Beside the tubular cast nephropathy, the CLL "accumulate" in the kidney even though with a lower frequency compared to MCN, in the light chains deposition disease (LCDD) and in the amyloidosis AL (AL). LCDD is characterized by a deposit of nodular amorphous materials PAS positive in the glomerulus and sometimes even in the tubulus. It usually presents itself as a chronic RI and a proteinuria causing nephrotic syndrome (NS). This quickly evolves into uraemia and its evolution can be lessened by the MM treatment. AL in course of MM also reeals with a chronic RI and NS. CLLs deposit in the typical fibrillar structure, on the vessel walls, in the glomerulus, in the mesangium and can be marked out with the Congo red colouring and the subsequent green birefringence through microscope with polarized light. Prognosis of AL is extremely severe and no benefit is given by the treatment of the hematological illness. It is therefore absolutely necessary to study the renal histology through biopsy when MM is grade B, that is, with serumal creatinine above 2 mg/dl as: MCN imposes the MM treatment programme in order to reduce the tubular excess of PBJ and to attempt to make RI reversible; MCN with tubular atrophy and interstitial fibrosis results in an unfavourable prognosis as it expresses a nephropathic irreversibility due to the loss nephrons. It will therefore necessary to start on a renal substitutional treatment programme. Renal damage in course of MM is not always tubular, rather an unexpected glomerular damage of LCDD or amyloidosis AL type can be found.(ABSTRACT TRUNCATED AT 400 WORDS).

23. Immunoturbidimetric assay for estimating free light chains of immunoglobulins in urine and serum.

J Clin Pathol 1991 Jun;44(6):466-71 (ISSN: 0021-9746)

Tillyer CR; Iqbal J; Raymond J; Gore M; McIlwain TJ

Department of Chemical Pathology, Royal Marsden Hospital, Sutton, Surrey.

An immunoturbidimetric assay for the assessment of free kappa and lambda light chains of immunoglobulins was developed using a commercial polyclonal antiserum with reactivity towards epitopes on the light chains, which are not expressed when they are bound to heavy chains. The assay, on a centrifugal analyser, is simple and rapid. The limit of detection is 5 mg/l of free light chain, with an assay range of 5-120 mg/l, intrabatch precisions from 1.5-6.4%, and interbatch precisions from 6.5-8.9%. The assay was only slightly less sensitive than colloidal gold staining of cellulose acetate electrophoreses for the detection of Bence-Jones protein in urine. For the serial monitoring of response to chemotherapy in patients with myeloma, the assay correlated well with serum paraprotein estimates obtained by densitometric scanning of Ponceau stained cellulose acetate electrophoreses, but not with serum beta-2 microglobulin measurements, even after correction for the effects of creatinine. These assays may prove to be of use for the monitoring of tumour response in the treatment of Bence-Jones myeloma.

24. Clinical applications of immunoglobulin free light chain analysis

Int J Clin Lab Res 1994;24(2):120-1 (ISSN: 0940-5437)

Pascali E

The potential of immunoglobulin-free light chain detection and measurement in biological fluids, and particularly in urine, has not yet been fully explored in clinical medicine, because of differences in sensitivity and lack of standardization of both quantitative assays and qualitative analysis. The ability to identify monoclonal free light chain, i.e., Bence-Jones protein, reliably in urine is of great clinical interest even when it occurs only in such small amounts that they frequently remain undetected by the routine methods used in many laboratories. The standardized adoption of both qualitative and quantitative procedures with comparable sensitivity is suggested to allow precise characterization of the clonal nature, and therefore the clinical role, of any free light chain excess.